

museu de ciències naturals de Barcelona Passeig Picasso s/n 08003 Barcelona T (+34) 93 256 22 22 museuciencies.cat

LPGyP RATES

*(A) *(B)

	Uncovered	Small	13	23
THIN SECTIONS		Medium	20	34.5
	Covered with protecting lacquer (Struers)	Small	16	27.5
		Medium	24	35
	Covered with coverslip	Small	18	31
		Medium	28	42
	With metallographic polishing	Small	25	43
	Carbonate staining		6	11
	Metallographic polishing only		11.5	20
	Coating removal (carbon, gold, etc)		2.5	5
	Metallographic specimens		30	55
POLISHING	Manual	Max 5 cm	10.5	18.5
		Max 15 cm	21.5	37.5
	With Al, Ce or Cr oxide		8.5	14.5
	With diamond. Add		9.5	17
	EpoFix epoxy resin (Struers) €/50 mL		15	27
	EpoThin low-viscosity epoxy resin (Buehler) €/50 mL		23	50
	EpoColor red epoxy resin (Buehler) €/50 mL		40	80
RESIN	Blue Pigment resin stain (Buehler) for porosity. Add €5/mL		20	35
	Polyester resin €/L		15	20
SAWING	Manual cutting with high-speed saw €/h		25	43
	Precision cutting with low-speed saw. By cut		4	7
	CRUSHING €/h (2)		28	49
	SIEVING €/h (2)		25	43
NEW	Smear slides		15	25
PROTOCOLS	Acetate peels		8	15
TECHNICIAN HOUR (3)			30	40

(A) Universities and public institutions

(B) Private companies

* Prices include VAT

(2) LPGyP self-service applications

(3) Rate only applicable in self-service applications if an LPGyP technician performs the work of the user in cleaning and/or finetuning equipment and facilities